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## Gepan® instill

**Please read the instructions carefully before use, as they contain important information about the use of Gepan® instill.**

Instructions for use

**Gepan® instill**

Sterile sodium chondroitin sulphate solution (0.2 %) for the temporary replacement of the glycosaminoglycan layer (GAG layer) of the bladder in cases of chronic cystitis, e. g. interstitial cystitis, overactive bladder, chronic recurrent urinary tract inflammation, radiation cystitis.

**Properties:**

The inner wall of the bladder is lined with an impermeable protective layer known as the GAG layer. This GAG layer is composed of a complex network of glycosaminoglycans and proteoglycans. The most prevalent constituent of the GAG layer is chondroitin sulphate. The main function of this layer is to prevent potassium ions in the urine, as well as micro-crystals, bacteria and proteases, from coming into contact with and irritating the underlying epithelium (outer layer of cells). Damage to the GAG layer interferes with this protective function. This allows irritating substances such as potassium ions, micro-crystals, bacteria or proteins to damage the bladder epithelium.

The bladder has some ability to repair its GAG layer by taking up endogenous urinary glycosaminoglycans into the ionically bonded network. In addition, exogenously supplied glycosaminoglycans can also be incorporated into the GAG layer as a temporary replacement in case of insufficient endogenous glycosaminoglycan production.

Gepan® instill was specifically developed for the temporary replacement of the GAG layer on the bladder epithelium and contains highly purified chondroitin sulphate, prepared using a special process. Gepan® instill is for instillation into the bladder only.

**Composition:**

40 ml solution contains 80 mg sodium chondroitin sulphate (0.2 %, 2 mg/ml). Gepan® instill contains no preservatives.

**Warnings:**

**Gepan® instill is for instillation into the bladder only.**

**Do not administer Gepan® instill systemically.**

**Gepan® instill is for single use only, as sterility cannot be guaranteed for a second use.**

**Discard unused portions.**

**Do not use Gepan® instill after the expiry date.**

**Also do not use if the disposable syringe is damaged or exhibits any signs of being unsuitable for use (e. g. clouding of the solution, liquid in the blister, cap not on the syringe).**

**Instructions for use:**

Instill the entire volume of Gepan® instill into the bladder by means of a catheter after the bladder has been completely emptied. For this, first remove the cap from the syringe and connect the syringe to an appropriate catheter using the **multi-catheter-compatible** connector.

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Gepan® instill should then remain in the bladder for as long as possible, at any rate for a minimum of 30 minutes. During the first four weeks of therapy, once-weekly instillation is recommended. Subsequently, monthly instillations should be conducted until symptoms have disappeared. Depending on the course of symptoms, some patients may require the initial phase of weekly instillations to be extended for a period of six or more weeks before changing to monthly instillations. Before use, warm the product to room temperature or, if necessary, body temperature. The presence of an air bubble in the syringe is due to technical reasons.

**Storage conditions:** Do not store above 25 °C. Do not freeze.

**Contraindications:**

Gepan® instill should not be used in cases of known hypersensitivity to the ingredients. Gepan® instill should not be used to treat pregnant women or children under 12 years of age, as there is no therapy experience in these patient groups.

**Interactions:**

Interactions with other medical devices for temporary replacement of the GAG layer are unknown. However, due to the lack of empirical experience, concurrent use with other GAG-layer-replacing medical devices should not take place. Interactions with medicinal products are unknown. However, due to theoretical considerations, concurrent treatment with medicinal products which act as urine acidifiers (e. g. methionine, ammonium chloride) should not take place.

**Side effects:** None known.


**Shelf life:** The expiry date is printed on the package and on the blister label. Do not use Gepan® instill after the expiry date.

Also do not use if the syringe is damaged.

**Pack sizes:** Prefilled syringe with 40 ml chondroitin sulphate solution.

**Manufacturer:** G. Pohl-Boskamp GmbH & Co. KG, Kieler Strasse 11, 25551 Hohenlockstedt, Germany

**Date of information:** September 2013

 25°C				
Temperature limit: Store between 2°C and 25°C!	See instruction for use!	Lot number	Do not inject!	For one-time application
				
Multi-catheter-tip	Expiry date	Sterile filtered under aseptic conditions	Manufacturer	

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**Produsent:** G. Pohl-Boskamp GmbH & Co. KG, Kieler Strasse 11, 25551 Hohenlockstedt, Tyskland