Test your health

PRIMA LAB SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND

support@primalabsa.ch
primalabsa.ch

EC REP Qarad EC-REP BV
Pas 257
2440 Geel - Belgium

800025_IFU_11_5.3 012024

REF 800025

Registered trademarks.

INSTRUCTIONS FOR USE GEBRAUCHSANWEISUNG NOTICE D'UTILISATION ISTRUZIONI PER L'USO

IRON FER TEST

Self-test for the detection of ferritin levels
in whole blood samples

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ENGLISH IRON AND ANAEMIA

Iron is an essential metal for our bodies and is crucial for transporting oxygen in the blood, for cell multiplication and to build the structure of tissues and organs. However, in excessively high levels it is toxic for the body. For this reason, every one of us has a system for taking up iron from the external environment (e.g. via a diet rich in iron-containing foods) and storing it in cells in a way that is not excessive (and therefore not toxic). Ferritin is the protein responsible for this storage function. The level of ferritin is an excellent indicator of the amount of iron available to the body. Low levels of this protein in the blood are an indication of depleted iron stores, a condition that precedes the development of anaemia. A decrease can be caused by anaemia, pregnancy, haemorrhages, alterations in iron uptake, tuberculosis.

TEST PRINCIPLE


Iron FER Test is an immunochromatographic assay which detects the protein Ferritin thanks to special monoclonal gold-conjugate antibodies embedded to test strip.

REFERENCES

1. Wick M, Pingerra W, Lehmann P. Iron metabolism: diagnosis and therapy of anemias, 5th ed, Vienna, New York: Springer Verlag, 2003; p. 151.
2. Worwood M. The laboratory assessment of iron status – an update. Clin Chim Acta 1997; 259: 3-23.
3. Kaltwasser JP, Werner E. Diagnosis and clinical evaluation of iron overload. Baillieres Clin Haematol 1989; 2; 363-89.
4. Baynes RD, Cook JD. Current issues in iron deficiency. Curr Opin Hematol 1996; 3:145-9.
5. Lee MH, Means RT Jr. Extremely elevated serum ferritin levels in a university hospital: associated diseases and clinical significance. Am J Med 1996; 98: 566-71.

CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

- * Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).
- * Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked.
- * The desiccant packet must not be used. Dispose of it with household waste without opening it.

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
EC REP Hualan Trianda Medical Instruments Co.,
Ltd., No.106 East Songjiang Road, Huaiyin,
Economic & Technological Development Zone,
223002 - Hualan City, Jiangsu - P. R. China

EC REP Phoenix Innovative Healthcare Manufacturers Pvt. Ltd.,
E1-209, Shi Mahape Road, Electronic Zone, MIDC,
TTC Industrial Area, Mahape, Navi Mumbai
400 710 MH 1 India

EC REP Shanghai International Holding Corp.
GmbH (Europe)
Eiffelstrasse 80,
20537 Hamburg - Germany

EC REP Advena Ltd., Tower Business Centre,
2nd Flr, Tower Street, Swatar,
BKR 4013, Malta

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
PRECAUTIONS

1. Read these instructions for use carefully before performing the test. The Test is reliable only if all the instructions are followed correctly.
2. Keep the Test out of the reach of children.
3. Do not use the Test after the expiry date or if the package has been damaged.
4. Follow the procedure exactly, using only the specified quantities of blood and diluent.
5. Store the Test components at a temperature between +4°C and +30°C. Do not freeze.
6. Use the test and lancing device once only.
7. The test is for external use only. DO NOT SWALLOW.
8. In vitro diagnostic device for individual use.
9. Not recommended for people who take anti-coagulant medications (blood thinners) or people suffering from haemophilia.
10. After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.

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







ENGLISH TEST PROCEDURE

- 1) Wash hands with soap and warm water, rinse with clean water and allow to dry. If this is not possible, use the provided gauze as an alternative.
2) Prepare the necessary material as follows: open the aluminium pouch, take out only the test cassette and throw away the desiccant bag.
3) Carefully rotate the protective cap of the sterile lancet 360° without pulling it.
4) Carefully massage the finger chosen for the puncture...
5) Holding the hand down, massage the finger until a large drop of blood forms.
6) Take the pipette without pressing the bulb.
7) Place the blood collected with the pipette into the well indicated on the cassette (S) by pressing the pipette bulb.
8) Unscrew the blue cap from the dropper vial.
9) Wait 5 minutes and read the results referring to the next section for result interpretation.

RESULTS INTERPRETATION

READ THE RESULT AFTER EXACTLY 5 MINUTES

The intensity of the line colours is not relevant for the purposes of interpretation of the Test's results.

NEGATIVE RESULT

Two coloured bands appear in the reading window by the T (test) and C (control) signs. The T band may be less intense (lighter) than the C line. This result means the levels of ferritin in the blood are above 30 ng/ml.

POSITIVE RESULTS

A coloured band appears only under the C (control) sign. This means the value of ferritin is below normal levels and medical advice is needed. Consult a physician.

NON VALID RESULT

No bands appear or there is a line only under the T (test) sign and not under the C (control) sign. In this case it is not possible to interpret the result of the test, which must be considered not valid. Repeat the test with a new blood sample.

F.A.Q. - QUESTIONS AND ANSWERS

HOW DOES THE IRON FER TEST WORK? Ferritin is a protein responsible for storing iron in cells. A positive result means that the concentration of ferritin is <= 30 ng/ml.

HOW DO I INTERPRET THE TEST IF THE COLOUR AND INTENSITY OF THE TEST AND CONTROL BANDS ARE DIFFERENT? The colour and the intensity of the bands are not important for the interpretation of the result.

WHEN CAN THE TEST BE USED? The Iron Fer Test can be performed in the presence of symptoms such as paleness, tiredness, frequent migraines and widespread pain, palpitations and occasionally increased heart rate, lack of muscle strength.

IS THE RESULT RELIABLE IF READ AFTER 10 MINUTES? No. The test must be read 5 minutes after the procedure is completed.

WHAT SHOULD I DO IF THE RESULT IS POSITIVE? If the result is positive, the concentration of ferritin in the blood is below normal levels and you should consult a physician, who will establish what to do next.

WHAT SHOULD I DO IF THE RESULT IS NEGATIVE? If the result is negative, the concentration of ferritin is at least 30 ng/ml. In some cases, samples with ferritin levels just below this value can yield a result of this type.

HOW ACCURATE IS THE IRON FER TEST? The Test is very accurate. Assessment reports show a concordance rate of 96.6% (CI 95%: 93.4-99.8%) with reference methods.

HOW ACCURATE IS THE IRON FER TEST? The Test is very accurate. Assessment reports show a concordance rate of 96.6% (CI 95%: 93.4-99.8%) with reference methods. Despite the reliability of the test, falsely positive or falsely negative results are possible.

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